

## SmartNPWT

Smart NPWT is wound care modality used in modern healing process. It is a treatment using sub-atmospheric pressure to increase blood flow, remove bacteria and increase growth of granulation tissue in the wound.

### SPECIAL FEATURES

- Light Weight, Portable Therapy Unit
- Provides Real – Time info on the therapy
- Vaccu-Sense Algorithm helps monitor & maintain target pressure
- Self Test mode to Identify & Resolve leaks
- Enhanced Fluid dynamics to minimise tubing blocks & alarms

### ORDERING INFORMATION

- Smart NPWT control unit LV 10010000
- Canister & Consumables kit LV 10020001

**MDE**  
MEDICAL DEVELOPMENT & ENGINEERING

## SMART NPWT NEGATIVE PRESSURE WOUND THERAPY

### Technical Specifications

<b>Negative pressure</b>	Vacuum from 20 - 150 mmHg adjustable in increments of 10 mmHg
<b>Voltages</b>	100-240 V ~ 50/60 Hz
<b>Sound emission</b>	<38 dB (A) (acc. to ISO 7779)
<b>Dimensions (WxHxD)</b>	9.8" x 10.2" x 4.2"
<b>Weight (Basic Unit )</b>	4.5 lbs.
<b>Ambient Conditions</b>	-30... +50°C, 95% RH

**MDE**

MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026  
Tel : +91-20-27121512, 30720117 Fax : +91-20-27121512  
Email : mdeindia@gmail.com Website : www.mdeindia.net

