



### INTELLIGENT<sup>®</sup>DERMATOME 2008

The Intelligent Dermatome 2008 consists of a precision engineered dermatome with wide selection of cutting widths and thicknesses combined with the performance of variable speed high torque motor unit for achieving very accurate and consistent cutting performance.

#### SPECIAL FEATURES AND BENEFITS

- ▶ Compact and light weight
- ▶ Wide selection of cutting depth from 0.20 to 0.38mm
- ▶ Powerful high torque motor unit in its class
- ▶ Dermatome heads, motor units and cable are sterilisable
- ▶ Widths of Graft from 25mm to 73mm
- ▶ Accurate and consistent cutting performance
- ▶ CE mark conformance
- ▶ Easily maintained and secured
- ▶ Micro computerized speed control for optimum torque output of motor

#### Technical Specifications - Dermatome Head

Cutting width	- 73 mm , 48 mm , 25 mm
Cutting thickness	- 0.20 to 0.38 in steps of 7 settings
Sterilization	- Autoclavable
RPM	- 5000-15000 rpm in steps of 100
Control	- Variable speed electronic control

## INTELLIGENT DERMATOME

### Micro Motor and Power Supply

Supply	- 230VAC
Classification according to MDD	- CE
Motor	- constant torque Brushless DC
Dimensions	- 225 (w) x 73 (h) x 199 (d) mm
Weight	- 1.5 Kgs

Note : Specifications are subject to change without notice

### Ordering information

Intelligent Dermatome 2008	LG 100 00000
Accessories	
Dermatome Head Block	LG 100 10000
Foot Switch	LG 200 20000
Motor with Cable	LG 100 30000
Blade	LG 100 40000
Support Plate 73 mm	LG 100 40001
Support Plate 48 mm	LG 100 40002
Support Plate 25 mm	LG 100 40003
Sterilization Container	LG 100 50000

## MDE

### MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026  
 Tel : +91-20-27121512, 30720117 Fax : +91-20-27121512  
 Email : mdeindia@gmail.com Website : www.mdeindia.com



## MDE

MEDICAL DEVELOPMENT & ENGINEERING

## INTELLIGENT DERMATOME

