

MI CRO MULTI TASK



- Electrical Drill System with power and high speed up to 30,000 rpm
- Console to allow visible display and setting of speed
- Foot control for varying the speed & foot button for reverse rotation
- System gives audible beeps/alerts while in reverse action
- Brushless DC motor for max torque & reliability
- Lightweight, flexible and autoclavable cables
- Quick release and lock system for tools and attachments
- Various blades options are available in Sagittal Saw
- Light weight and non-bulky motor
- Powerful motor with power 120W & Torque 25mN-m to reduce vibration
- No lubrication or seal required to run the motor

MI CRO MULTITASK

Micro Multitask System



Hand Piece
Part No : LM 10020000



Control Unit
Part No : LM 10000000

Common Attachments



Drilling attachment with
mini Jacob's Chuck
Part No : LM 10030001



Keyless Dedicated Drill 1.5 mm
Part No : LM 10030002



Keyless Dedicated Drill 2.2 mm
Part No : LM 10030003



Keyless Dedicated Drill 3.2 mm
Part No : LM 10030004



Micro Sagittal Saw
Part No : LM 10050001



Reciprocating Saw
Part No : LM 10060001

MDE

MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026
Tel : +91-20-27121512, 30720117 Fax. : +91-20-27121512
Email : mdeindia@gmail.com Website : www.mdeindia.net

