

**POWERDRIVE
N - SHAVER**

MDE
MEDICAL DEVELOPMENT & ENGINEERING

**POWERDRIVE
N - SHAVER**



Shaver Handpiece
Part No : LS 10020000



Full Radius Resector
Part No : LS 10050001



Agressive Full Radius
Part No : LS 10050002



Meniscus Cutter
Part No : LS 10050003



Round Burr
Part No : LS 10050004



Tapered Burr
Part No : LS 10050005



Oval Burr
Part No : LS 10050006



Corded Powerdrive
Part No : LC 10020000



Drilling with Jacob's Chuck
Part No : LB 10020001



Reaming AO Synthes Coupling
Part No : LB 10020002



Sagital Saw
Part No : LB 10020003



Mini Sagital Saw
Part No : LB 10020004



Wire Driver
Part No : LB 10020005

MDE

MEDICAL DEVELOPMENT & ENGINEERING

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026 India.
Tel : +91-20-27121512, 30720117 Fax. : +91-20-27121512
Email : mdeindia@gmail.com



POWERDRIVE N - SHAVER



POWERDRIVE N-SHAVER Features

- Introducing " Power - N - Shaver " in one console
- Suitably designed to provide adequate torque & speed to all arthroscopic , drilling and saw applications
- Loaded with TFT touchscreen, provides surgeons feedback on speed & operating mode.
- Can attach shaver Handpiece or Modular Handpiece to one unit.
- Modular Handpiece can accommodate drilling, saw, reamer and wire driver attachments like any operating drill.
- Shaver attachment can accommodate — shaver blades to suit applications.




POWERDRIVE N - SHAVER



Footswitch
Part No : LS10030000



Control Unit
Part No : LS 10010000



Shaver Handpiece
Part No : LS 10020000



Corded Powerdrive
Part No : LC 10020000

