

Battery Operated

POWERDRIVE  
MULTITASK



Weight 1.4 kg without attachment

POWERDRIVE MULTITASK  
Motor Features

- Cannulated shaft.
- Integrated gear for optimum transmission of torque.
- Specially designed hand-crafted brushless servo motor for maximum torque in its class.
- Autoclavable, Integrated 3ø servo drive ensures maximum torque and excellent reliability.
- Ergonomic design ensures easy manoeuvrability during all procedures.

Battery Operated

POWERDRIVE  
MULTITASK

Battery Operated Systems

Handpiece for Drill & Reamer



Hand Piece  
Part No : LB 10010000

Battery Casing



Battery Charging Unit  
Part No : LB 10040000

Battery Status Indicator  
Charging time 90 min.



Asceptic Battery Packs  
Part No : LB 10050000  
Life - 500 charge cycles



Sterile Transfer Shields  
Part No : LB 10060000

Dedicated Saw Handpiece



Hand Piece  
Part No : LB 10010001

Features

- Heavy duty saw for orthopaedic procedure
- Multiple position rotatable head
- 14500 cpm, 3.5°
- Various sizes of blades for osteotomy THR, TKR

## Battery Operated

### Common Attachments



Drilling with Jacob's Chuck  
Part No : LB 10020001  
Cannulation 4.0mm  
0.4 to 6.7 mm Drill bit



Reaming AO Synthes type  
Part No : LB 10020002  
for Intramedullary Reaming  
0-270 rpm

### Added Attachment



Reciprocating Saw  
Part No : LB 10020008

\* These products are not part of standard supply.  
Special request to be made for quotation.

# MDE

**MEDICAL DEVELOPMENT & ENGINEERING**

34, Electronic Sadan 1, MIDC, Bhosari, Pune 411 026

Tel : +91-20-27121512, 30720117 Fax. : +91-20-27121512

Email : mdeindia@gmail.com

## POWERDRIVE MULTITASK



Mini Sagital Saw  
Part No : LB 10020004  
14500 cpm, 3.5°



Wire Driver  
Part No : LB 10020005  
for wire ø 1-4mm



Acetabular Reaming  
Part No : LB10020006  
0-270 rpm



Quick Coupling  
Part No : LB10020007

Accessories like Sagital Saw Blades, Reamers,  
Drill Bits, could be requested extra.



# MDE

MEDICAL DEVELOPMENT & ENGINEERING



## POWERDRIVE MULTITASK

